

AfriDX

COVID-19 diagnostics for Africa

Autoinduction Protein Expression of Engineered Boon2-BST

Cambridge Analytical Biotechnology Group
Department of Chemical Engineering and Biotechnology
University of Cambridge



E D C T P



This project was funded by EDCTP and the European Union's Horizon 2020 research and innovation program under grant number RIA2020EF-2918

Preliminary note to user:

- The following protocol was designed using BL21 (DE3) as the host *E. coli* strain
- All the steps of the “protein expression protocol” section should be performed under sterile conditions until step 11 (with the exception of the incubation steps)

Materials:

- Sodium phosphate dibasic (Na_2HPO_4)
- Potassium phosphate monobasic (KH_2PO_4)
- Sodium chloride (NaCl)
- Tryptone
- Yeast extract
- Glycerol
- Glucose
- Lactose
- Water
- LB Broth
- Bunsen Burner or Biosafety cabinet (for environment sterility)
- 2 autoclavable 1 litre bottle
- 3 autoclavable 250 mL bottles
- Conical Flask
- Centrifuge Tubes
- Pipettes
- Kanamycin
- Glycerol stock of BL21 (DE3) transformed with the pET24a – Boon2-BST plasmid

Prepare autoinduction media:

Prepare the following for autoclaving.

Solution 1 (2x Media): Weigh out and add the following to a 1 litre autoclavable bottle:

- 6 g of Na_2HPO_4
- 3 g of KH_2PO_4
- 5 g NaCl
- 20 g of Tryptone
- 5 g Yeast Extract
- 500 mL of water



Solution 2 (50% vol/vol glycerol): Pour the following into a 250 mL autoclavable bottle:

- 50 mL of glycerol
- 50 mL of water

Solution 3 (10% weight/vol glucose): Add/Pour the following into a 250 mL autoclavable bottle:

- 5 g of glucose
- 100 mL of water

Solution 4 (5% weight/vol lactose): Add/Pour the following into a 250 mL autoclavable bottle:

- 10 g of lactose
- 200 mL of water

Solution 5: Pour at least 500 mL of water into a 1 L autoclavable bottle.

Autoclave all 5 solutions.

Reconstitute autoinduction media: Under sterile conditions add the following to the sterilised 2x media (solution 1):

- 12 mL of sterile 50% glycerol
- 5 mL of sterile 10% glucose
- 40 mL of sterile 5% lactose
- Top up to 1 L using sterile water

Protein Expression Protocol of Engineered Boon2-BST:

The protein expression process starts with an overnight culture.

Step 1: Using 70% ethanol, sterilise the benchtop/biosafety cabinet to be used.

Step 2: Using 70% ethanol, sterilise and clean any pipettes that are to be used.

Step 3: Add 5 mL of LB Broth to a sterile centrifuge tube.

Step 4: Add kanamycin to the LB Broth for a final concentration of 50 µg/mL. (I.e. if the stock concentration of the kanamycin is 100 mg/mL, add 2.5 µL of the stock to 5 mL of LB Broth.)

$$\frac{50 \frac{\mu\text{g}}{\text{mL}} \times 5 \text{ mL}}{100 \frac{\text{mg}}{\text{mL}}} = 2.5 \mu\text{L}$$



Step 5: Inoculate the prepared LB Broth with 1 μL of the Boon2-BST bacterial glycerol stock.

Step 6: Incubate the culture overnight (12 to 16 hours) at 37 °C and 225 rpm.

The following steps are performed the next morning.

Step 7: Add the desired culture volume worth of autoinduction media to a sterile conical flask (i.e. 50 mL, 100 mL, 200 mL, etc). The maximum volume is dictated by the volume of the flask. Standard protocols utilise an air-to-volume ratio of 10:1. This means that if a 2 L flask is being used, the max desired volume of LB Broth should be 200 mL.

Step 8: Add kanamycin to the autoinduction media for a final concentration of 50 $\mu\text{g}/\text{mL}$. (i.e. if the stock concentration of the kanamycin is 100 mg/mL and the desired culture volume is 200 mL, add 100 μL of the stock to the LB Broth.)

$$\frac{50 \frac{\mu\text{g}}{\text{mL}} \times 200 \text{ mL}}{100 \frac{\text{mg}}{\text{mL}}} = 100 \mu\text{L}$$

Step 9: Inoculate the prepared culture flask with the overnight culture at a ratio of 100:1. (i.e. if the desired culture volume is 200 mL, add 2 mL of the overnight culture.)

Step 10: Incubate the large culture at 37 °C and 225 rpm for 10 hours.

Step 11: Aliquot the culture equally into centrifuge tubes.

Step 12: Centrifuge the culture at 3750 rpm for 30 minutes.

Step 13: Decant the supernatant. The pellets can be stored at -20 °C long term or 4 °C short term.